

# KERATOMETER

SD-125 - Speedway, Delhi (INDIA)




since 1934

- L **Coincidence Focusing System - Positive measurement upto 0.02 mm accuracy**
- | **15 X Wide Angle Piece**
- | **Measures both meridians without changing optical system**
- | **Simple vertical adjustment to fit patient**
- | **Two - way adjustable head & chin rests**
- | **Dual Eye Level Sighting System, facilitates horizontal alignments**
- | **Positive fixation permits rapid measurements of central corneal area.**
- | **Precision objectives assures durability**
- | **Achromatic Lenses**
  
- | **Measuring Range**
  - o **Range of Curvature: 9.4 ~ 6.4 mm (in 0.05 mm steps)**
  - o **Diopter: 36 ~ 52 (in 0.25 steps)**
  
- | **Standard Accessories**
  - o **Dust Cover**
  - o **Spare Bulb**
  
- | **Optional Accessories**
  - o **Instrument Table Manual / Motorized**



Specifications of equipment/s are subject to change without prior notification

ISO 9001:2008;  
ISO 13485:2003 




since 1934

**Speedway**  
**SURGICAL CO.**

*Precision is our vision*

E-mail: [sales@speedwaydelhi.com](mailto:sales@speedwaydelhi.com) Website: [www.speedwaydelhi.com](http://www.speedwaydelhi.com)

 +91-11-2386 6379 / 4563 0593 / 2722 7545 / 2741 7545 / 4237 2101