

LENSOMETER

SD-175 - Speedway, Delhi (INDIA)



since 1934



Type	:	Ocular, External Reading Type
Target	:	Pinhole Corona & Cross Target
Vertex Power Range	:	(0.25 D Steps Upto 10D & 0.50D Steps Beyond 10D)
Cylindrical Axis	:	0° - 180° (In 1° Increments)
Prism Diapter Scale:	:	5 Prism Diopters (In 1d Increments)
Lens Size	:	Upto 90 mm
Eye Piece Adjustment	:	Possible
Tilting	:	Possible
Power Supply	:	3v Dc (2 Pencil Cells)
Weight	:	4.5 Kgs
Size	:	20 x 12 x 22 (L x B x H)

Specifications of equipment/s are subject to change without prior notification

ISO 9001:2008; ISO 13485:2003 



since 1934

Speedway
III SURGICAL CO.

Precision is our vision

E-mail: sales@speedwaydelhi.com Website: www.speedwaydelhi.com

 +91-11-2386 6379 / 4563 0593 / 2722 7545 / 2741 7545 / 4237 2101