

OPERATION TABLE

SDOT-2000 - Speedway, Delhi (INDIA)



SPEEDWAY

since 1934

Key Features :

Vertical Height Adjustment

Enables the Ophthalmologists to perform from A convenient, comfortable standing or sitting Positions

Optimally Positioned Head Rest

Manually controlled up & down option allows The surgeons for easy access to the patient While offering the patient comfortable head Positioning.

Tredelenberg

Built in Trendelenberg (Head Low) & Reverse Trendelenberg (Head up) Positions ensure Safety During emergencies

Misc. Features

Ample Leg Room Sufficient space for equipment Pessals Easy to clean & maintain Convenient Stability during operations

Ergonomic Operation Table with comprehensive accessories rang

→ IV Pole

→ Built in Oxygen Inlet Tube
→ Slide Tray


Technical Specifications

Height Range min. / max.	:	640 / 870 mm
Lifting speed	:	8 mm / sec
Total length	:	1880 mm
Total width	:	690 mm
Base size	:	825 / 520 mm
Input Definition	:	230 V 50/60 Hz. 110 V 50/60 Hz.
Weight Capacity	:	135 Kgs.
Gross Weight	:	115 Kgs.

Tilting Possible: 18°
Manual Head Adjustment
Available with Safety Belt and
two 5 Amp sockets
Height Adjustment with foot padle



Specifications of equipment/s are subject to change without prior notification

ISO 9001:2008; 
ISO 13485:2003



SPEEDWAY
since 1934

Speedway
III SURGICAL CO.
Precision is our vision

E-mail: sales@speedwaydelhi.com Website: www.speedwaydelhi.com

+91-11-2386 6379 / 4563 0593 / 2722 7545 / 2741 7545 / 4237 2101