

REFRACTION UNIT, Manual

SD-2001 - Speedway, Delhi (INDIA)



Features :

One fully upholstered elegant ophthalmic chair

1. With full manual recline facilities

One Stand Console:

1. With Illuminating soft light for examination
2. With controls for ophthalmic chair
3. For "Recline " and "Up & Down"
4. For Up & Down Movement pedal id also provided

Specification :

Base dimensions

Floor space required

- | | | |
|--------------------------|---|-------------------------------|
| 1. Height | : | 7.6" (Appox.) |
| 2. Length | : | After reclining 9.6" (Appox.) |
| 3. Width | : | 4" (Appox.) |
| 4. Input Voltage & Power | : | 110/12V |
| 5. Output | : | 3V,6V, 12V |
| 6. Weight | : | 240kg. (Appox.) |

Optional Accessroies :

1. Slit lamp
2. Binocular Indirect ophthalmoscope
3. Keratometer or Autorefractometer
4. Direct Ophthalmoscope
5. Chair Projector
6. Trail Lens Set
7. Sight Tester (Phoropter)



Specifications of equipment/s are subject to change without prior notification

ISO 9001:2008;
ISO 13485:2003



SPEEDWAY
since 1934

Speedway

III SURGICAL CO.

Precision is our vision

E-mail: sales@speedwaydelhi.com Website: www.speedwaydelhi.com

+91-11-2386 6379 / 4563 0593 / 2722 7545 / 2741 7545 / 4237 2101