

TRIAL LENS SET

S-8321 - Speedway, Delhi (INDIA)



since 1934



S-8321GS:01 Moorfield Trial Case Reduced Aperture Illuminated (Golden & Silver)
S-8321GS:02 Moorfield Trial Case Reduced Aperture Non Illuminated (Golden & Silver)

Specifications:

- Trial Case, with reduced aperture lenses having color Metal rims for quick identification
- 32 pairs each spherical lenses CX & CC 0.12D to 20.0D & cylindrical 18 pairs each of CX & CC 0.25 to 6.0.
- Seven Accessories - one each Maddox rod, Slit, Blank, Pinhole, Red, Green & Plain Prisms-1/2, 1, 2, 3, 4, 5, 6, 8, 10, 12
- Spherical 20D (every .25D up to 4D, every 0.5D up to 6D, every 1D to 14D and every 2D to 20) Cylinder - 6D (every 0.25D of up to 3.5D, every 0.5D upto 6D)

S-8321RB:01 Moorfield Trial Case Reduced Aperture Illuminated (Red & Black)
S-8321RB:02 Moorfield Trial Case Reduced Aperture Non Illuminated (Red & Black)

Specifications:

- Trial Case, with reduced aperture lenses having color Metal rims for quick identification
- 32 pairs each spherical lenses CX & CC 0.12D to 20.0D & cylindrical 18 pairs each of CX & CC 0.25 to 6.0.
- Seven Accessories - one each Maddox rod, Slit, Blank, Pinhole, Red, Green & Plain Prisms-1/2, 1, 2, 3, 4, 5, 6, 8, 10, 12
- Spherical 20D (every .25D up to 4D, every 0.5D up to 6D, every 1D to 14D and every 2D to 20) Cylinder - 6D (every 0.25D of up to 3.5D, every 0.5D upto 6D)



S-8321F:01 Moorfield Trial Case Full Aperture Non Illuminated (Small)

Specifications:

- Trial Case, with full aperture lenses
- 25 pairs each spherical lenses CX & CC 0.25D to 20.0D & cylindrical
- No Accessories
- No Prism
- Spherical 20D (every .25D up to 3D, every 0.5D up to 5D, every 1D upto 12D and every 2D upto 20D) Cylinder - 6D (every 0.25 of up to 3D, every 1D upto 6D)

S-8321F:02 Moorfield Trial Case Full Aperture Non Illuminated (Large)

Specifications:

- Trial Case, with full aperture lenses
- 40 pairs each spherical lenses CX & CC 0.25D to 20.0D & cylindrical
- Accessories Cross Cylinder 2, Ocluder 1, Red 1, Green 1, Plain 2, Slit 2, Pinhole 2, Frosted 1, Bagolini 1 (Total of 13 Accessories)
- Prism 0.5, 1, 2, 3, 4, 5, 6, 8, 9, 10 (Total of 11 Prisms)
- Spherical 20D (every .25D up to 6D, every 0.5D up to 10D, every 1D upto 16D and every 2D upto 20D) Cylinder - 6D (every 0.25 of up to 4D, every 0.5 upto 6D)



Specifications of equipment/s are subject to change without prior notification

ISO 9001:2008;
ISO 13485:2003



since 1934

Speedway
SURGICAL CO.
Precision is our vision

E-mail: sales@speedwaydelhi.com Website: www.speedwaydelhi.com

+91-11-2386 6379 / 4563 0593 / 2722 7545 / 2741 7545 / 4237 2101